IZI Medical Products
 5 Easter Court, Suite J
 Owings Mills, MD 21117
 USA
 Tel: +1 410 594 9403

IZI Medical

Rx only

DPA-PAK (Disposable Passive Array – Pedicle Access Kit MDT)

Pedicle access needle with cannula (1)

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Navigated Pedicle Access Kit (EN)

Set Components

- Pedicle access needle with cannula (1)
- Removeable Handle (2)
- Extra cannula (3)

Indications for Use

The Disposable Passive Array - Navigated Pedicle Access Kit is intended to be used in image guided surgery with a Medtronic StealthStation® System. It allows an accurate placement of pedicle screws. The device is sterile and designed for single use.

Warnings

- All components in the Pedicle Access Kit are pre-sterilized, single-use components. Do not re-sterilize or reuse. Attempts to re-process or re-sterilize the device may be ineffective and may compromise its structural integrity. Any re-use creates a risk of contamination which could result in patient injury, illness, or death.
- Inspect the Set components for any damages, shelf life expiration, deterioration and abuse prior to procedure. If found, do not use the device.
- The Product is MR UNSAFE.
- The product must only be used by or under the guidance of a skilled, duly qualified physician trained in the use of a Medtronic StealthStation system manufactured by Medtronic Navigation Inc.
- This operating manual is not a reference to procedural technique.
- Do not modify the device. Any modification may compromise device safety and efficacy.
- Surgeons must determine the appropriate use of the Pedicle Access Kit.
- During navigation, visually confirm navigational accuracy frequently by localizing on known anatomical points, including accuracy checkpoints, and comparing the position of the instrument in the image with its physical location in surgical space.
- The Pedicle Access Kit components are designed for use only when fully assembled. Make sure that the handle, needle, and cannula are assembled together and that the cannula is securely locked in place.
- If you disconnect and reconnect the Pedicle Access Kit components (for example, to change cannulas or for another approach), re-calibrate the device before navigating.
- Use care when handling the needle because it has a sharp tip that may cause injury.

Cautions

- Federal Law (U.S.A) restricts this device to sale by or on the order of a physician.
- Refer to the Medtronic application instructions for use for general navigation instructions, warnings, and precautions.

Handling and Application Information

1. Verify that the Navigated Pedicle Access Needle is enabled in the application software.
2. Calibrate the needle by placing the tip of the needle in the reference frame divot and pressing the system foot switch or holding the needle in the divot until the Medtronic StealthStation system responds.
3. Use the application software to navigate the needle into position for the procedure using standard surgical technique.
4. Properly dispose of the kit components after use. Do not attempt to re-sterilize or re-use these items.

Assistance

If you have any questions about the Pedicle Access Kit, contact IZI Medical Products +1 410 594 9403 for support.

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