Thermalfix™

Radiopaque bone cement

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Directions for use

First EC marking date applied: 2006

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner) or on the physician's prescription. Before using IZI Medical products, the operator must thoroughly familiarise himself with the safety regulations figuring in the directions for use, as well as the information relating to each product (description, operating technique, leaflets, etc.). The corresponding information is available from IZI Medical. The surgeon should also be informed of the residual risk associated with the product he considers using.

GENERAL INFORMATION

Implantation of IZI Medical products should only be performed by qualified operators with in-depth knowledge and perfect command of the specific operating techniques of IZI Medical products. Operative techniques may be acquired from the manufacturer or the distributors. The surgeon is responsible for any complications or harmful consequences which may arise due to an incorrect indication or operating technique, improper use of the material and failure to observe the safety regulations provided in the directions for use. Neither the IZI Medical manufacturer nor the authorized IZI Medical representative can be held responsible for these complications. This product is a fast-setting acrylic cement for vertebroplasty or kyphoplasty. It is available in the form of a sterile liquid ampoule and a sterile powder pouch (liquid monomer and powder polymer).

COMPOSITION

Powder (27,2g) :	
Polymethylmethacrylate	49.5 %
Benzoyl peroxide	0.5 %
Zirconium dioxide	45.0 %
Hydroxyapatite	5.0 %
Liquid (9,2g) :	
Methyl Methacrylate	99 %
N-N dimethyl-p-toluidine	1 %
Hydroquinone	20 ppm

INJECTION SET

S5 injection kit (sold separately)
- One bevelled trocar 12,5 cm in length and 11G (sold separately)

INDICATIONS

The Thermalfix™ is indicated for the fixation of pathological fractures of the vertebral body using vertebroplasty or kyphoplasty procedures. Painful vertebral compression fractures of the vertebral body may result from osteoporosis, benign lesions (hemangioma), or malignant lesions (metastatic cancers, myeloma).

CONTRAINDICATIONS

The use of the Thermalfix™ osseous cement is contraindicated in patients having one of the following problems:
- Cases of active, suspected or incompletely treated infection

- Coagulation disorders, or severe cardiopulmonary disease
- Spinal stenosis (>20% by retroplused fragments)
- Compromise of the vertebral fracture due to posterior involvement
- Patient clearly improving on more conservative medical treatment
- Prophylaxis in metastatic or osteoporotic patients with no evidence of acute fracture
- Non-pathological, acute traumatic fractures of the vertebra
- Vertebral plana (collapse>90%)
- Compromise of the vertebral body or the walls of the pedicles
- Unstable vertebral fractures due to posterior involvement
- Hypersensitivity to one of the constituents of the product.

ADVERSE EVENTS

Serious adverse events, some with fatal outcome, associated with the use of acrylic bone cements for vertebroplasty or kyphoplasty include: myocardial infarction cement, cardiac arrest, cerebrovascular accident, decrease in blood pressure, pulmonary embolism, sudden death, short-term cardiac conduction disorders and cardiac embolism. Although the majority of these adverse events present early within the post-operative period, there have been some reports of diagnoses beyond a year or more after the procedure.

Other reported adverse events for acrylic bone cements intended for vertebroplasty or kyphoplasty include:

- Leakage of the bone cement beyond the site of its intended application with introduction into the vascular system resulting in embolism of the lung and/or heart or other clinical sequelae.
- pneumonia, intercostal neuralgia, pneumothorax, fracture of a pedicle;
- fracture of ribs in patients suffering from diffuse osteopenia, in particular during thoracic vertebroplasty procedures, due to the high pressure exerted downwards while the needle is being inserted;
- collapse of a vertebra adjacent to the injected vertebra due to an osteoporotic disease,
- leak of cement into the intervertebral disks;
- leak of cement into the vascular system
- leak of cement into soft tissues
- leak of cement with compression of the spinal cord possibly resulting in paralysis or loss of sensitivity;

Interactions with other agents: none known to date.

INSTRUCTIONS FOR USE

1- Setting up the trocars

Once the track has been determined by the operator, the trocars are set in place under radioscopic guidance, using a digitizing table with a hoop.

2- Preparation of Thermalfix™ cement

Powder:

Open the peel-off bag and take out the sterile powder pouch

Carefully open the pouch and pour all the whole of the powder into the injection kit mixer.

Liquid:

Open the ampoule blister pack and take out the sterile ampoule.

Do not break the ampoule over the mixer (risk of glass splinters).

Pour the whole of the liquid onto the powder.

Mix for 45 seconds

If the S5 injection kit is used, follow the directions for use mentioned in the instructions for this kit.

3- Filling the injector

- Transfer the cement into the cement injector.

The low viscosity of Thermalfix™ cement facilitates the operation.

Note : when mixing with a bowl and a spatula, mix for 45 seconds. Cover the bowl to prevent evaporation of monomer. Wait two minutes and stir the cement slowly for 10

seconds. The cement is then ready to be used.

4- Setting up and injection

Screw the injection device onto the trocar pre-positioned in situ. The injection of the cement should be performed under continuous radiological guidance. In order to prevent vascular migration, the cement should be inserted during its injection phase.

Injection is to be stopped when the operator estimates that the vertebral filling is satisfactory, or when a risk of leak of cement becomes apparent.

Once the body of the vertebra has been filled, insert the mandrel in the trocar used for the injection to prevent cement residues from setting in soft tissues when the trocar is withdrawn. At a 20°C temperature of the operating room and of the equipment, the different phases are as follows:

With a Gauge 11 trocar:

- Mixing : 0'45»
- Filling : 0'45»- 3' 00» *
- Injection time : 3' 00» - 12' 30» *
- Setting : 12' 30» - 24' 30» **

*times obtained with the S5 injection system. They may vary if different systems are used.

**test based on the ISO 5833 standard

These times are more than sufficient for the operator to conduct continuous radiological controls and space out the fillings, thereby avoiding any unwanted migration of the cement.

!! WARNING !!

- The length of phases of application depends on room temperature and that of components, but also on the degree of hygrometry in the operating block. A high temperature reduces hardening time. A low temperature extends this time.

- For controlled and optimal use of the Thermalfix™ cement, the doses are to be stored at the intended operating room temperature (recommended 20°C) for a minimum of 24 hours before use.

- Read instructions carefully prior to use.

- The operator should have specific training and experience to be thoroughly familiar with the properties, handling characteristics, application of the product and percutaneous cement delivery.

- The manufacturer does not recommend a surgical technique: it is the practitioner's responsibility to determine the appropriate use of Thermalfix™ bone cement and specific technique for each patient.

- Follow cement handling, mixing and preparation instructions carefully.

-It is strictly forbidden to re-sterilize products. Single use only. Sterile only if package is unopened and undamaged.

- Hypotensive reaction have occurred, between 10 to 165 seconds. They have lasted from 30 seconds to 5 minutes Some patients have progressed to cardiac arrest. For this reason, patients should be monitored carefully for any change in blood pressure during or immediately after the application of bone cement.

- Methyl methacrylate may cause hypersensitivity among high-risk persons, which can result in an anaphylactic reaction.

- Due to insufficient data, the safety and efficacy of this cement has not been established with regard to pregnant women and children.

The addition of HA has not been shown to provide any improved clinical results.

- We do not recommend using these products on patients that do not suffer from a pathological condition, such as primary or secondary osteoporosis or a tumour. This could impair the ability of the patient to recover using conservative treatment methods.

- Established surgical principles and techniques must be strictly observed. Deep wound infection is a serious postoperative complication that may require complete removal of the bone cement. Deep wound infection may be latent and may not become apparent until several years after the operation.

- Caution should be exercised to prevent excessive exposure to the concentrated fumes of the monomer, which may irritate the respiratory tract, eyes and even the liver.

- Always check the condition of the liquid before carrying out the procedure. Do not use the liquid component if it shows any sign of thickening or premature polymerisation. These conditions indicate that the product has not been stored correctly.

- Do not allow the liquid component to come into contact with rubber or latex gloves. The liquid component is a powerful lipid solvent. Should contact occur, the gloves might dissolve, causing tissue damage. Wearing two pairs of gloves may reduce the possibility of hypersensitivity reactions.

- Do not allow personnel wearing contact lenses to be near or involved in mixing of the bone cement.

- Use appropriate imaging techniques to verify that the needle is correctly positioned, that no damage has been caused to surrounding structures, and that the injected bone cement has been correctly located. Use an imaging technique, such as fluoroscopy, to assess the capacity of the vertebra to contain the injected bone cement.

- Avoid over-pressurizing the bone cement, as this may cause the bone cement to leak beyond the site of its intended application. Cement leakage may cause tissue damage, and nerve or circulatory problems. Excessive pressure can also be applied to the injection device, which could cause it to break.

- Leaks may also occur during the injection if the needle is in a vein or if there are undetected micro-fractures.

- If bone cement is detected outside the vertebral body or in the circulatory system during the procedure, stop the injection immediately.

- Inadequate fixation or unanticipated postoperative events may affect the cement / bone interface and lead to micromotion between the cement and bone surface. A layer of fibrous tissue may then develop between the cement and the bone. Long-term regular supervision if therefore recommended for all patients.

- The final polymerisation stage occurs in situ and is an exothermic reaction with considerable liberation of heat. According to ISO 5833 standard, the temperature can reach 95°C. Maintain the patient's position until the polymerization process has been completed, so as to obtain proper fixation. An additional 1 to 2 hours or more may be necessary, depending on the patient's medical condition and the operator.

- The long term effects of the bone cement in the spine have not been established.

- Cement leakage may cause tissue damage, nerve or circulatory problems, and other serious adverse events.

PRECAUTIONS FOR USE

The surgeon must be trained in Thermalfix™ use technique and must follow it scrupulously. We therefore advise you observe the preparation times recommended by the manufacturer.

Adverse effects may occur whenever you fail to comply with the technique guide.

The application phases of Thermalfix™ cement have been determined using an injection system recommended by IZI Medical. They may vary depending on the injection system used.

- It is also recommended that the mixing times given should be complied with in order to prevent the polymerization from progressing to such an extent that the cement is no longer fluid enough to allow easy transfer into the injection system, and proper filling of the bone cavities.

During the application of Thermalfix™, radiological control is essential so that the operator can follow the progress of the filling and stop the procedure if the slightest leakage of cement is detected.

A thorough preoperative check-up of the patient must be carried out prior to the operation.

- Make sure the operating room is properly ventilated to eliminate monomer fumes as much as possible.

The monomer is a volatile and flammable liquid.

Ignition of monomer fumes caused by the use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported.

- The insertion of a foreign body into the tissues increases the normal risk of infection associated with surgery during the postoperative period.

This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

STERILIZATION

- The liquid in the ampoule is sterilized by ultra-filtration, and the ampoule blister is sterilized using ethylene oxide. The powder is inside a double sterile pouch. This double pouch is sterilized by 25 kGy gamma rays.

- Prior to using, carefully check the protective packaging in order to ensure that it has not been damaged in a way that could affect its sterility.

- When taking the product out of its packaging, be sure to follow the asepsis rules.

- The cement is delivered sterile and ready for use in the operating room.

- Re-sterilization of the product is prohibited. Single use only

- Do not use after the use-by date.

- Disposable. Sterile only if package is unopened and undamaged.

IMPORTANT INFORMATION TO PHYSICIANS

- Percutaneous vertebroplasty or kyphoplasty procedures should only be performed in medical settings in which emergency decompressive surgery is available.

- Adverse reactions affecting the cardiovascular system have been attributed to acrylic cement. Recent data indicate that the monomer undergoes rapid hydrolysis into methacrylic acid, and that a significant fraction of the circulating methacrylate is present in the form of free acid rather than methyl ester. Correlation between changes in circulating concentrations of methyl methacrylate/methacrylic acid and changes in blood pressure has not been established.

- Check for any change in blood pressure during and immediately after application of the osseous cement.

- The doctor is responsible for any complication or harmful consequences, which may result from an erroneous indication or operating technique, from inappropriate use of the material, or non-observation of the safety instructions that figure in the directions for use.

- Additives (such as antibiotics) are not to be mixed with the bone cement, as this will alter cement properties.

INFORMATION TO THE PATIENT

The patient should be informed by the doctor of the potential consequences of the factors mentioned in the following paragraphs: contra-indications and adverse events, that is, those liable to hinder the success of the operation, as well as possible complications which may arise. The patient should also be informed of the measures to be taken to diminish the possible consequences of these factors.

PACKAGING STORAGE

Designation	Powder (g)	Liquid (g)
Thermalfix™	27.2g	9.2g

The cement is to be stored in its original, unopened packaging, in a dry, clean place, away from light and at a temperature of 25°C max.


Keep away from any ignition source.

RECOMMENDATIONS FOR DISPOSAL

-Allow the cement to harden before disposal with other medical waste. Comply with local regulations in force relating to medical waste for safe handling and disposal of the cement.

-Regarding separate disposal of the liquid or of the powder, comply with local regulations in force relating to handling and disposal of the cement.

**Influence Of Temperature On Handling Times :
(+/- 1 min)**

* times obtained with the TEKNIMED S5 injection system
 ** test based on the ISO 5833 standard

- Catalogue number
- Batch code
- Manufacturer
- Do not reuse
- Upper limit of temperature
- Flammable product
- Consult instructions for use
- Caution
- Date of manufacture YYYY-MM-DD
- Do not resterilize
- Use by date YYYY-MMMM-DD
- Keep away from sunlight
- Do not use if package is damaged
- Product compound, using alternative way of sterilization
- Sterilized using irradiation

R_X ONLY CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Manufactured in France by:
TEKNIMED S.A.S
 8, rue du Corps Franc Pomiès
 65500 VIC EN BIGORRE - FRANCE
 Tél. (33)5 62 96 88 38
 Fax (33)5 62 96 28 72
 www.teknimed.com

Distributed by:
IZI Medical Products
 5 Easter Court, Suite J
 Owings Mills, MD 21117
 USA
 Tele: +1 410 594 9403
 Fax: +1 410 594 0540
 www.izimed.com

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