



Vertefix[®] MV

Bone Cement Infused with Insite™ Tracking Beads

Used for fixation of pathological fractures of the vertebral body using vertebroplasty and kyphoplasty procedures.



- Two components (powder & liquid) PMMA bone cement
 - Cement components are easy to mix, no measurement is required
 - Injection time: 18 minutes. Yield: 18 cc of PMMA
- Supplied sterile in box. Intended for one-time use.

Indications of use are for the fixation of vertebral compression fractures resulting from:

- Malignant lesions (Metastatic cancers, myeloma)
- Benign lesions (Hemangioma)
- Osteoporosis and osteopenia



Scan To Learn More

Optimal Visibility and Versatility.

ORDER INFORMATION

Vertefix[®] MV

PRODUCT NAME	MIXED VOLUME	CEMENT RADIOPACITY	VISCOSITY	GLOBAL PRODUCT NUMBER
Vertefix MV	18cc	30%	Medium	G90006

Kit Configurations

PRODUCT NAME	DESCRIPTION	GLOBAL PRODUCT NUMBER
VCF Kit	DuroJect Injector, Minimix, and Vertefix MV	G80009
10G Bone Filler Kit	Osteoforce Injector, 10G Bone Fillers, Syringes, Minimix, & Vertefix MV	G80010
11G Bone Filler Kit	Osteoforce Injector, 11G Bone Fillers, Syringes, Minimix, & Vertefix MV	G80011

Not all parts numbers shown may be approved for sale in all regulatory jurisdictions. Parts and configurations are subject to change without notification. Please consult with your Account Representative for details.